


## Jay - Zip (Paediatric)



- Dual layered contoured foam base
- Reticulated foam in the outer cover promotes air movement
- Soft and stretchable outer cover
- All models include: 2 way stretch cover

### Warranty Information

- 2 Years Cushion
- 6 Months Cover

 Contract approved

### ndis

**Item Code:**

05\_181006190\_0105\_1\_2

**Pressure Care Cushion –**

**Composite:**05\_429\_0103\_1\_2

Click [here](#) for more information

# JAY<sup>®</sup>

### Compliance

No applicable compliance standard

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC585184	JAY ZIP - 10 x 10in	2.75 in	75 kg	Medium Risk	Foam
PTC585193	JAY ZIP - 12 x 12in	2.75 in	75 kg	Medium Risk	Foam
PTC585185	JAY ZIP - 12 x 16in	2.75 in	75 kg	Medium Risk	Foam
PTC585192	JAY ZIP - 14 x 14in	2.75 in	75 kg	Medium Risk	Foam
PTC585190	JAY ZIP - 14 x 16in	2.75 in	75 kg	Medium Risk	Foam
PTC585191	JAY ZIP - 15 x 15in	2.75 in	75 kg	Medium Risk	Foam