

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	WAS	NOW	QTY
BEA007481	Alarm - Cura - Bed Monitor Plug Pack 9V DC (replacement)	\$ 46.20	\$ 13.86	1
BEA007520	Alarm - Monitor - Mono - To Suit Nurse Call Systems	\$ 184.00	\$ 55.20	4
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	5
BEA024126	Bed Extender - Invacare - 100mm - for Invacare All CS Beds - Adjustable Width	\$ 550.00	\$ 165.00	5
BEA028100	Overbed/Chair Table - Aspire - Wide C Base - Thermoform Recessed Split Tilt Top - Cream	\$ 349.00	\$ 174.50	1
BEA030050	Bed Mover/Transporter - Aidacare FL250 Bed (Set of 2)	\$ 809.00	\$ 404.50	1
BEBO45100	Bed - Long Double - 4 Function - 400kg	\$ 8,719.00	\$ 2,615.70	1
BEBO45850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	1
BEM046725	Mattress - Aspire Lifecomfort Premier Care - Single - 180kg	\$ 869.00	\$ 434.50	1
BTC062900	Shower Commode Seat - K Care - H/Made - Open 550mm - Blue	\$ 579.00	\$ 173.70	1
BTC077000	Pan Carrier - for K Care Junior 400mm Shower Commode - S/Steel	\$ 293.00	\$ 87.90	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	1
LAC353105	Linen Bag - Half Size (OHS) - White	\$ 18.70	\$ 5.61	4
LAC353106	Linen Bag - Half Size (OHS) - Blue	\$ 18.70	\$ 5.61	5
LAC353107	Linen Bag - Half Size (OHS) - Green	\$ 18.70	\$ 5.61	5
LAC353108	Linen Bag - Half Size (OHS) - Red	\$ 18.70	\$ 5.61	5
LAC353109	Linen Bag - Half Size (OHS) - Yellow	\$ 18.70	\$ 5.61	5
LAC353110	Linen Bag - Half Size (OHS) - Purple	\$ 18.70	\$ 5.61	5
LIC360870	Ultra Control Unit - incl, hand control - Fits Likorall ES and rails	\$ 1,629.00	\$ 488.70	1
LIC360900	UltraTwistTM Wide - incl, hand control. Max 500 kg (1100 lbs).	\$ 3,829.00	\$ 1,148.70	1
LSS371615	Lifter Sling - Oxford - Quickfit Poly Padded - Large	\$ 459.00	\$ 229.50	3
LSS371698	Lifter Sling - Oxford - Quickfit Deluxe Poly - Small	\$ 217.00	\$ 108.50	1
LSS373400	Sling General Purpose Fab Small	\$ 327.00	\$ 163.50	1
LSS380500	Lifter Sling - Gen Purp Mesh - Kerry - Small LR406NE1	\$ 297.00	\$ 148.50	2
LSS380700	Lifter Sling - Gen Purp Mesh - Kerry - Large LR406NE3	\$ 297.00	\$ 148.50	2
LSS381800	Lifter Sling - Gen Purp Mesh & Head Support - Kerry - Small LR407NE1	\$ 349.00	\$ 174.50	1
LSS390918	Lifter Sling - Aspire - LGE - Deluxe General Purpose w/Head Supp - CLIPS - Mesh - 300kg	\$ 423.00	\$ 211.50	1
LSS390940	Lifter Sling - Aspire - SML - Deluxe Hygiene/Access w/Head Supp - CLIPS - Poly - 200kg	\$ 480.00	\$ 240.00	1
LSS390960	Lifter Sling - Aspire - SML - Low Friction General Purpose w/Head Supp - CLIPS - Parasilk - 200kg	\$ 539.00	\$ 269.50	3
LSS394120	Sling - SupportVest w/ padding Mod 91 - Liko - Sml - w/2 Loop Clips Fits Sabina II 350 SlingBar.	\$ 709.00	\$ 212.70	1
LSS394180	Sling - ComfortVest Mod 95 - Liko - Sml - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 779.00	\$ 233.70	1
LSS394200	Sling - ComfortVest Mod 95 - Liko - Lge - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 829.00	\$ 248.70	1
PTC585149	Cushion - JAY - Lite 17x17	\$ 482.00	\$ 241.00	1
PTC589500	Cushion - ROHO Quadro Select LP 20x16.75in (11x9cell)	\$ 1,119.00	\$ 559.50	1
PTC593900	Cushion - ROHO Cont Sel 20x18.5in (11x10cell)	\$ 1,119.00	\$ 559.50	1
PTC594946	Cushion - ROHO Hybrid Elite 1R 22x18in Sensor Ready	\$ 1,239.00	\$ 619.50	1
PTC594951SR	Cushion - ROHO Hybrid Elite 1R 20x20in w/HD cover Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594963	Cushion - ROHO Hybrid Elite 2R 16x20in	\$ 1,079.00	\$ 539.50	1
PTC608040	Cushion - Promedicare Domino - 16 x 16in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC608056	Cushion - Promedicare Domino - 19 x 18in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTM636391	Mattress - Foam - Invacare - Softform Premier - King Single - 1970x1070x150mm	\$ 1,679.00	\$ 839.50	1
WAT907201	Avalo AX Articulating Arm with Pole Mount (needs 1x A&E Mount 12149K)	\$ 447.00	\$ 134.10	1

AGED CARE
**WAREHOUSE
CLEARANCE**

ALL PRODUCTS
BETWEEN
50-70% OFF

**HURRY,
WHILST
STOCK
LASTS!**

Item Number	Item Description	WAS	NOW	QTY
-------------	------------------	-----	-----	-----