

Seating

Standing

Sleeping

Therapy

JENX[®]
Positioning for life



Multistander 2 by Jenx

A highly versatile 3 in 1 standing system; offering upright, prone and supine positioning for children from approx. 4 - 13 years.



Always Supportive...

Jenx Limited, Wardsend Road, Sheffield S6 1RQ, United Kingdom
tel: +44 (0) 114 285 3376
email: export@jenx.com
www.jenx.com

JENX[®]
Positioning for life



Multistander 2 by Jenx

This fully adjustable and versatile 3 in 1 standing frame is suitable for children aged from approx. 4 - 13 years; providing a high level of postural support in upright, prone or supine.

An unrivalled amount of adjustment for height and width and a maximum user weight of 60kg, make the Multistander 2 the perfect choice for growing children in multi-use environments. Smooth and easy angle adjustment takes the Multistander from horizontal right through to upright, facilitating easy transfers in supine and unlimited scope for setting the precise standing angle required for every child.

The wide variety of head, trunk, leg and foot support options enable you to exactly customise the Multistander 2 to any individual user, whilst the deliberately open design ensures you can always see precisely how the user is positioned each and every time. The capacity for outstanding support, coupled with the exceptional growth and the versatile range of accessories, mean the Multistander 2 will easily meet the changing needs of a growing child for many years.

INNOVATIVE
FEATURE

De-rotational Thoracic and Pelvic Straps with Lateral Support Pads

Straps fasten in the centre and adjust independently from either side, enabling a mid-line position to be achieved and maintained. Buckles and hook & loop fastener mean that once adjusted correctly the straps can be fixed, making it difficult for siblings or classmates to tamper with! **Support Pads** available in small, medium and high rise.

INNOVATIVE
FEATURE

Head Support

A range of different head support options are available to accommodate different support requirements, from mild support offered from the **Moulded PU Head Support** to the **Multigrip™ Headrest** for those who require much greater support.

INNOVATIVE
FEATURE

Prone/Supine Adjustable Tray

Adjustable for height, angle and depth, whilst unique inbuilt infills prevent those awkward gaps that arms get stuck in when standing in supine. Supplied with mounting slots for **Elbow Blocks** and a raised contrasting tray edging which helps provide both a visual barrier for children and a useful anchor point for attaching toys or other sensory aids.

INNOVATIVE
FEATURE

Multigrip™ Body Support

Can be used as a shoulder protractor in supine, or to provide additional lumbar support between the trunk and pelvic boards.

INNOVATIVE
FEATURE

Angle Adjustable Footplate

Provides customised foot positioning by allowing easy adjustments of dorsi/plantar flexion for those who cannot maintain plantigrade.





**INNOVATIVE
FEATURE**

**Independently
Adjustable Thoracic &
Pelvic Support Boards**

Allowing the individual proportions of each child to be catered for, whilst also giving scope to avoid pressure on delicate areas such as feeding buttons.

**INNOVATIVE
FEATURE**

Angle Adjustment

The Multistander 2 offers the option of either a powered or manual base to adjust the tilt-in-space angle. The powered base removes any manual handling concerns and is controlled by a hand-held controller. The manual base has a robust winder mechanism which provides a cost effective alternative. When used in supine, the frame can be adjusted to horizontal to also aid with transfers.

**INNOVATIVE
FEATURE**

Practical Modern Look

The white frame provides a fresh modern look, whilst making it easy to spot spills of body fluid/feed and clean before they become hazardous.

Always Encouraging...

Standing is a great way to improve your child's outlook on the world and their physical development.

Improves blood pressure & circulation



Aids respiration & speech

Helps digestion, bowel & bladder emptying



Helps relieve pressure created during sitting

Encourages bone & muscle development



Supports pelvis, thighs & feet so children can join in everyday activity

Gives a new perspective and ways to participate



Provides opportunities to increase social skills



Children also gain huge psychological advantages from standing, as well as those important physical and development benefits.

INNOVATIVE FEATURE

PU Support Pads

Easy to clean and impervious to fluids! The boards and **Support Pads** on the Multistander 2 are made from polyurethane that has an integral antimicrobial agent within it. This ensures Multistander 2 is quick and easy to keep clean and hygienic, which aids with re-issue for stores or for meeting infection control for multiple users.



Always Configurable Accessories

Multistander 2 by Jenx is a modular product, allowing you to choose specific accessories to meet the exact requirements of each child.

Multistander 2 Frame

The Multistander frame includes:

- Manual or powered base options
- Base with 4 lockable swivel castors
- Knee block mounting brackets*
- Single angle adjustable footplate

*NB: **Knee Blocks** must be specified separately.



Support Boards

PU chest & pelvic **Support Boards** in a range of 5 colours.



Colour options available.



Multigrip™ Headrest (supine standing only)

The original and innovative Jenx **Multigrip™ Headrest** is available in 3 sizes and offers adjustment for height, depth and angle.



Oval or Moulded PU Pad Headrest (supine standing only)

Headrests for children who require minimal head support. Comes mounted on a flexible fixing bracket, giving adjustment for height depth and angle.



Colour options available.



Support Pads

PU **Support Pads** can be used at the hips or chest and are adjustable for height, width and angle and can be off-set to brace a scoliosis or accommodate and support unique positioning needs.

Support Pads come complete with appropriate size de-rotational strap as standard.



Colour options available.



Flexible Thoracic Support Pads

Available in 2 sizes, these super flexible support pads (with straps) are designed to provide outstanding positioning with comfort in mind. Independently adjustable for height and width to ensure maximum trunk control.



Cover colour options available.



Multigrip™ Body Support

Can be used as a shoulder protractor in Supine, or to provide additional lumbar support between the trunk and pelvic boards



Cover colour options available.



Waistcoat Harness

Provides additional trunk and shoulder stability for children with low tone. An innovative zip system which zips from the bottom up allows for easy-access to gastronomy feed tubes. Available in 3 sizes with a wipe-clean outer surface fabric.



Tray

Tray for use in prone or supine. Easy to attach and adjust for height, angle and depth, complete with a removable tray infill for use in prone standing.



Chest Prompt

Small foam sternal prompt, used to encourage extension of the thorax and help improve head control. Requires **Support Board Covers** for use.



Support Board Covers

Removable covers which wrap around the centre of the chest and pelvic **Support Boards** to add a splash of colour and for use with **Chest Prompt**.



Colour options available.



Tray Elbow Blocks

Used to encourage arms into a midline position and/or to block excessive retraction.



Knee Blocks

Innovative and flexible PU **Knee Blocks** available in 2 sizes, adjustable for height, width, angle and depth. Complete with padded, de-rotational knee straps.



Colour options available.



Sandals

Foot positioning available in 2 sizes to accommodate different size feet, even when wearing splints. Jenx **Sandals** are interchangeable across the entire range of Jenx seating and standing products.



Ankle Snugs

An optional foot accessory to the standard **Sandals**, the **Ankle Snugs** offer a unique magnetic buckle fastening mechanism, which can be operated single-handedly.



Sandal Raising Blocks

Available for every size of Jenx **Sandals**, the raising block kits give the **Sandals** that extra bit of height where needed! Use to accommodate leg length discrepancy or to raise the height of the footplate to maximise the available growth.



NB: One kit needed per individual sandal.

Multistander 2 Colour Options



Blue



Pink



Green



Red



Grey

Visit jenx.com for more information.



**INNOVATIVE
FEATURE**

Knee Blocks

Available in 2 sizes and provided with specially designed straps, the **Knee Blocks** are independently adjustable for height, width, depth and angle, allowing you to achieve the optimum position for each child.

Multistander 2 by Jenx Technical Specifications

IMPORTANT INFORMATION:

Please be aware, some minimum and maximum measurements are achieved by the use of additional accessory components.

	Measurement	Units	Multistander 2 MUS02
Age Range	approx.	yrs	4yrs - 13yrs
User Height	min - max	mm	1100 - 1550
User Weight	max	kg	60
Chest Width	min - max	mm	200 - 320
Hip Width	min - max	mm	200 - 320
Height (sandal plate to top of upper support board)	min - max	mm	750 - 1200
Height (sandal plate to top of head support: Supine only)	min - max	mm	910 - 1350 (no shoulder support) 910 - 1500 (with shoulder support)
Height (floor to top of upper support pad)	min - max	mm	830 - 1280
Kneeblock Width (centre to centre)	min - max	mm	135 - 210
Footplate Angle	dorsi - plantar	degrees	15° - 15°
Tilt-In-Space Angle	horizontal - upright	degrees	0° - 85° (Powered base) 10° - 85° (Manual base)
Base Footprint	width x length	mm	675 x 1070

Find Your Distributor

To locate one of our carefully selected distributors near you, visit:
jenx.com/where-to-buy



@JenxLtd



@JenxLtd



@JenxLtd



Jenx Limited, Wardsend Road, Sheffield S6 1RQ, United Kingdom
tel: +44 (0) 114 285 3376
email: export@jenx.com
www.jenx.com

JENX[®] 
Positioning for life